NaCl 0,9% B. Braun

BG	Инструкции за употреба	2
CZ	Návod k použití	3
DE	Gebrauchsanweisung	4
DK	Brugsanvisning	5
EE	Kasutusjuhend	6
ES	Instrucciones de uso	7
FI	Käyttöohje	8
FR	Mode d'emploi	9
GB	Information for use	10
GR	Πληροφορίες σχετικά με τη χρήση	11
HR	Upute za uporabu	12
HU	Használati utasítás	13
IT	Istruzioni per l'uso	14
LT	Vartojimo instrukcija	15
LV	Lietošanas instrukcija	16
MK	Упатство за користење	17
MY	Maklumat penggunaan	18
NL	Gebruiksaanwijzing	19
NO	Brukerveiledning	20
PL	Instrukcja zastosowania	21
PT	Instruções de uso	22
RO	Informații pentru utilizare	23
SE	Bruksanvisning	24
SI	Navodilo za uporabo	25
SK	Informácia pre používateľa	26

schwarz

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Lätus - - -

Medical Device_Irrigation Solution
- Spuellösungen____190_GIF
190/12609811/0120
GIF (A5-Booklet, 25-spr.)
Production site: Crissier

Font size: 8,0 pt.
7,0 pt. (Address, p.18)

V-0100

B | BRAUN

- 1 -

NaCl 0,9% B. Braun

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0116

NaCl 0,9% B. Braun

Sodium Chloride 0,9% for Irrigation
Irrigation solution,
sterile and endotoxin-free

Composition

1000 ml solution contain:

Sodium chloride 9.00 g
in Water for Injection

Electrolytes:	mmol/l
Na ⁺	154
Cl ⁻	154

Indications

For irrigation and cleansing during surgical intervention.
For irrigation of wounds and burns.
For moistening wound tamponades, cloths and dressings.
For rinsing and cleansing of bladder catheters.
For irrigation and cleansing in stoma care.
For rinsing and cleansing of instruments and equipment.
For intra- and postoperative irrigation (endoscopic interventions without HF current, after TUR interventions).
For filling respiratory air humidifiers (the use of container equal or higher than 100 ml is recommended).
For mechanical irrigation of the eyes (the use of 30 ml container is recommended).
For nasal irrigation (the use of 30 ml container is recommended).

Dosage

The amount of fluid depends on the clinical condition.

Overdosage

Overdosage during wound irrigation and intra-operative application can lead to absorption of irrigation fluid, that can lead to hypervolaemia in the form of hypotonic hyperhydration (hyponatraemia).
Headache, nausea, restlessness and disorientation are symptoms of hyperhydration. In severe cases coma-like states can occur.

Countermeasures in cases of overdosage

Discontinuation of irrigation, correction of the water and electrolyte balance depending on the particular clinical situation.

Warnings and precautions

Do not use for infusion.

Only use the solution if the seal is undamaged and the solution is clear.

Opened containers must not be saved, quantities left over must be discarded. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Store out of reach of children. When there is a risk of absorption of larger amounts of the solution, as may be the case when using it in large body cavities or large wound surfaces, it should not be used in patients in states of hyperhydration. Under those conditions it should be used with caution in patients with

- Hypertonic dehydration,
- Hypokalaemia,
- Hypernatraemia,
- Hyperchloraemia,
- Disorders where restriction of sodium intake is indicated, such as cardiac insufficiency, generalised oedema, pulmonary oedema, hypertension, eclampsia, severe renal insufficiency.

NaCl 0,9% B. Braun must not be used for intrathecal irrigation procedures.

Adverse side effects

None are known.

Expiry date

The solution must not be used after the expiry date given on the label.

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B. Braun Melsungen AG